

SHIELD^{me}TM OXYGEN REGULATOR

CE



SHIELD^{me}TM OXYGEN REGULATOR controls the rate of continuous flow of oxygen, while an oxygen conserving device delivers a pulse of oxygen on demand to a patient as they breathe.

SPECIFICATIONS:

- Brass & Plastic
- Flow rate: 0 – 10/15lpm
- Pressure Gauge Protected Float Type
- Input Pressure: 150bar;
- Output Pressure: 2bar-3bar;
- Control Pressure of Safety Valve: 3.5bar± 0.5bar;
- Flowmeter Control Range: 0-10LPM/15LPM;
- Material: Main brass + aluminum alloy flow meter;
- Connecting thread: Female CGA540-type;
- With 200ml humidifier(reusable);
- Accessory: Nasal cannula attached with every carton.
- CE Approved

**We reserve the right to change specifications of our product without prior notice.*

info@shieldmeglobal.com

www.shieldmeglobal.com