

# SHIELD<sup>me</sup>™

## PROTECTING GOGGLE / EYE MASK

SME-MGG05

CE



**SHIELD<sup>me</sup>™ PROTECTING GOGGLE / EYE MASK** are developed for anti-virus, bacterial blocking rate over 95%, very suitable for blocking virus infection and achieve the preventing effect. The best choice to protect you and your family.

This high-quality goggle is waterproof and impact-resistant, not only can blocks saliva droplets, splash, and dust but also effectively isolate viruses. The frame is made of lightweight flexible PC. Feature adjustable, elastic straps to comfortably fit children and adults, the wraparound construction provides a perfect fit and seal. The protective goggles include 4 white indirect vents that allow maximum airflow. It has superior clear polycarbonate lens, full-view visual effect, and effective UV blocking. There is enough space inside our stylish goggles, which will not affect you wearing glasses.

### PRODUCT PERFORMANCE FEATURES

- ✔ This product is made of polymer material, light weight and strong, and can be used for many scenes, such as laboratory, hospital, outdoor and so on, can effectively prevent dust impact, liquid spatter or splash.
- ✔ Air-conditioning valves are provided on both sides to ensure safety and Suitable, Airtight during usage and ready to breathe through the valve. Lens with anti-fog coating, can effectively prevent the lens from fogging.
- ✔ The lenses are made of high-quality PC material with high clarity and ensure that they do not cause human body normal vision any influence, can wear to work with glasses at the same time.

- ✓ Ergonomic nose care design, no increased nose burden, more comfortable.
- ✓ Adjustable valve design on both sides, can be breathable at any time to prevent the influence of fog surface.
- ✓ Made of polymer materials to effectively prevent foreign matter impact and liquid spatter.
- ✓ The lens is made of high quality PC material, high light transmission and high clarity.
- ✓ Suitable for correct vision, and easy to adjust the head band, suitable for all kinds of head size.

## APPLICATIONS:



ANTI-IMPACT



ANTI-SPLATTER



LIGHTWEIGHT



HIGH LIGHT  
PENETRATION

*\*We reserve the right to change specifications of our product without prior notice.*